

## MEDEA Positive Control for STJ501750 peptide (STJ504451)

### GENERAL INFORMATION

|              |              |
|--------------|--------------|
| Product Type | Biomolecules |
| Applications | WB           |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                                                                          |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Dilution Range      | WB: 1:500-1:750                                                                                                                                                                                                                                          |
| Formulation         | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                                                                                                          |

### TARGET INFORMATION

|             |                                                                                               |
|-------------|-----------------------------------------------------------------------------------------------|
| Specificity | This is positive control is recommended for use in combination with MEDEA antibody STJ501750. |
|-------------|-----------------------------------------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|