

Anti-ATG4D antibody (265-474) (STJ110152)

STJ110152

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	84971
Gene Symbol	ATG4D
Uniprot ID	ATG4D_HUMAN
Immunogen	
Immunogen Region	265-474
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 265-474 of human ATG4D (NP_116274.3). DVTRLVVVVSQDCTVYKADV ARLVARPDPTAEWKSVILV PVRLGGETLNPVYVPCVKEL LRCELCLGIMGGKPRHSLYF IGYQDDFLLYLDPHYCQPTV DVSQADFPLESFHCTSPRKM AFAKMDPSCVTGFGYAGDRKE FETLCSELTRVLSSSSATER YPMFTLAEGHAQDHSLLDLC SQLAQPTLRLPRTGRLLRAK RPSEDFVFL
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081