

Anti-HTT antibody (550-650) (STJ110209)

STJ110209

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3064
Gene Symbol	HTT
Uniprot ID	HD_HUMAN
Immunogen	
Immunogen Region	550-650
Specificity	A synthetic peptide corresponding to a sequence within amino acids 550-650 of human Huntingtin (NP_002102.4).
Immunogen Sequence	LNDGTQASSPISDSSQTTTE GPDSAVTPSDSSEIVLDGTD NQYLGQLGIGQPQDEDEEATG ILPDEASEAFRNSSMALQQA HLLKNMSHCRQPSDSSVDKF V