

Anti-PIP5K1A antibody (231-500) (STJ110250)

STJ110250

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8394
Gene Symbol	PIP5K1A
Uniprot ID	PI51A_HUMAN
Immunogen	
Immunogen Region	231-500
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 231-500 of human PIP5K1A (NP_001129109.1).
Immunogen Sequence	MNNLLPRSVKMHKYDLKGS TYKRRASQKEREKPLPTFKD LDFLDIPDGLFLDADMYNA LCKTLQRDCLVLQSFKIMDY SLLMSIHNIHDAQREPLSSE TQYSVDTRRPAPQKALYSTA MESIQGEARRGGTMMETDDHM GGIPARNSKGERLLLYIGII DILQSYRFVKLEHSWKALV HDGDTVSVHRPGFYAERFQR FMCNTVFKKIPCVHLGRPDV LPQTTPLEEISEGSPIDPD



Western blot analysis of extracts of HeLa cells, using PIP5K1A antibody (STJ110250) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Western blot analysis of extracts of HeLa cells, using PIP5K1A antibody (STJ110250) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ3000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081