

## Anti-RNF125 antibody (1-232) (STJ110289)

STJ110289

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	54941
<b>Gene Symbol</b>	RNF125
<b>Uniprot ID</b>	RN125_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	1-232
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-232 of human RNF125 (NP_060301.2).
<b>Immunogen Sequence</b>	MGSVLSTDSGKSAPASATAR ALERRRDPPELPVTSFDCAVC LEVLHQPVRTRCGHVFCRSC IATSLKNNKWTCPCYCRAYLP SEGVDPATDVAKRMKSEYKNC AECDTLVCLSEMRAHIRTQ KYIDKYGPLQEEETAARCV CPFCQRELYEDSLDHCITH HRSERRPVFCPLCRIPDEN PSSFSGSLIRHLQVSHTLFY DDFIDFNIIEALIRRVLDLDR SLLEYVNHNSNTT