

Anti-USP26 antibody (664-913) (STJ110306)

STJ110306

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	83844
Gene Symbol	USP26
Uniprot ID	UBP26_HUMAN
Immunogen	
Immunogen Region	664-913
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 664-913 of human USP26 (NP_114113.1). TSLCQFHKAGGKPASSPGTP LSKVDFQTPENPKRKKYVK TSKFVAFDRIINPTKDLYED KNIRIPERFQKVSEQTQQCD
Immunogen Sequence	GMRICEQAPQQALPQSFPKP GTQGHTKNLLRPTKLNQKS NRNSLLALGSNKNPRNKDIL DKIKSKAKETKRNDKGDHT YRLISVSHLGLTKSGHYI CDAYDFEKQIWFTYDDMRVL GIQEAQMQRDRCTGYIFFY MHNEIFEMLKRENAQLN



Western blot analysis of extracts of various cell lines, using USP26 antibody (STJ110306) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Western blot analysis of various lysates using USP26 Rabbit polyclonal antibody (STJ110306) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081