

Anti-LAMB3 antibody (18-230) (STJ110426)

STJ110426

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:100
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3914
Gene Symbol	LAMB3
Uniprot ID	LAMB3_HUMAN
Immunogen	
Immunogen Region	18-230
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 18-230 of human LAMB3 (NP_000219.2). QQACSRGACYPPVGDLLVGR TRFLRASSTCGLTKPETYCT QYGEWQMKCKCDSRQPHNY YSHRVENVASSSGPMRWWQS
Immunogen Sequence	QNDVNPVSLQLDLDRRFQLQ EVMMFEQGMPAGMLIERSS DFGKTWRVYQYLAADCTSTF PRVRQGRPQSWQDVRCQSLP QRPNARLNGGKVQLNMDLV SGIPATQSQKIQEVGEITNL RVNFTRLAPVPQR