

Anti-DYM antibody (1-250) (STJ110521)

STJ110521

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	54808
Gene Symbol	DYM
Uniprot ID	DYM_HUMAN
Immunogen	
Immunogen Region	1-250
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DYM (NP_060123.3).
Immunogen Sequence	MGSNSSRIGDLPKNEYLKKL SGTESISENDPFWNQLLSFS FPAPTSSSELKLLLEATISV CRSLVENNPRTGNLALIKV FLSRTKELKLSAECQNHIFI WQTHNALFIICLLKVFICQ MSEEELQLHFTYEKSPGNY SSDSEDLLEELLCLMLQIT DIPLLDITYEISVEAISTMV VFLSCQLFHKEVLRQSISHK YLMRGPCLPYTSKLVKTLTY NFIRQEKPPPPGAHVFPQQ