

Anti-GTF2H2C_2 antibody (1-260) (STJ110596)

STJ110596

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	728340/730394
Gene Symbol	GTF2H2C.GTF2H2C_2
Uniprot ID	T2H2L_HUMAN
Immunogen	
Immunogen Region	1-260
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human GTF2H2C_2 (NP_001035955.1).
Immunogen Sequence	MDEEPERTKRWEGGYERTWE ILKEDESGSLKATIEDILFK AKRKRVEHHGQVRLGMMRH LYVVVDGSRTMEDQDLKPNR LTCTLKLLEYFVEEYFDQNP ISQIGIIVTKSKRAEKLTEL SGNPRKHITSLKEAVDMTCH GEPSLYNSLSMAMQTLKHMP GHTSREVLIIFFSLTTC DPS NIYDLIKTLKAAKIRVSVIG LSAEVRVCTVLARETGGTYH VILDESHYKELLTHLSPP



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081