

Anti-YTHDF3 antibody (1-275) (STJ110693)

STJ110693

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:100-1:500 IF/ICC:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your speci
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	253943
Gene Symbol	YTHDF3
Uniprot ID	YTHD3_HUMAN
Immunogen	
Immunogen Region	1-275
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human YTHDF3 (NP_689971.4).
Immunogen Sequence	MSATSVDRPKGGGNKVSQ NGSIHQKDAVNDDDFEYLS SQTNQSNYSPPMSDPYMPSY YAPSIGFPYSLGEAAWSTAG DQPMPLYLTYYGQMSNGEHY IPDGVFSQPGALGNTPPFLG QHGFNFFPGNADFSTWGTSG SQGQSTQSSAYSSSYGYPPS SLGRAITDGGAGFGNDTLK VPGISSIEQGMTGLKIGDDL TAAVTKTVGTALSSSGMTSI ATNSVPPVSSAAPKPTSWA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081