

Anti-EXO5 antibody (1-373) (STJ110805)

STJ110805

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64789
Gene Symbol	EXO5
Uniprot ID	EXO5_HUMAN
Immunogen	
Immunogen Region	1-373
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-373 of human EXO5 (NP_073611.1).
Immunogen Sequence	MAETREEETVSAEASGFSDL SDSEFLEFLDLEDAQESKAL VNMPGPSSSESLGKDDKPISL QNWKRGLDILSPMERFHLKY LYVTDLATQNWCELQATYQK ELPGLAPEKAAVLDTGASI HLARELELHDLVTPVTTKE DAWAIKFLNILLIPTLQSE GHIREFPVFGEGEVLLVGV IDELHYTAKGELELAELKTR RRPMLPLEAQKKKDCFQVSL YKYIFDAMVQGVTPASLI



Western blot analysis of extracts of various cell lines, using EXO5 antibody (STJ110805) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081