

Anti-NAA25 antibody (703-972) (STJ110814)
STJ110814

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	80018
Gene Symbol	NAA25
Uniprot ID	NAA25_HUMAN
Immunogen	
Immunogen Region	703-972
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 703-972 of human NAA25 (NP_079229.2).
Immunogen Sequence	LNHPVEPKNSEKTAENGVSS RIDILRLLLQQLLEATLETGK RFIEKDIQYPFLGPVPTRMG GFFNSGCSQCQISSFYLVND IYELDTSGLEDTMEIQERIE NSFKSLLDQLKDVFSKCKGD LLEVKGDNLKTHTPTLLENLV FVETISVILWSSYCESVL RPYKLNLQKKKKKKKETSII MPPVFTSFQDYVTGLQTLIS NVVDHIKGLETHIALKLEE LILEDTSLSPEERKFSKTV



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081