

Anti-Phospho-MKNK1-T197/202 antibody (STJ110872)

STJ110872

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8569
Gene Symbol	MKNK1
Uniprot ID	MKNK1_HUMAN
Immunogen	ITTPELTT
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around T197 & T202 of rat MKNK1 (NP_001037732.1).
Immunogen Sequence	ITTPELTT