

Anti-MAPKAPK3 antibody (1-240) (STJ11100068)

STJ11100068

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7867
Gene Symbol	MAPKAPK3
Uniprot ID	MAPK3_HUMAN
Immunogen	
Immunogen Region	1-240
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human MAPKAPK3 (NP_004626.1).
Immunogen Sequence	MDGETAEEQGGPVPPVAPG GPGLGGAPGGRREPKEYAVT DDYQLSKQVLGLGVNGKVLE CFHRRTGQKCALKLLYDSPK ARQEVDPHWHQASGGPHIVCI LDVYENMHGKRCLLIIMEC MEGGELFSRIQERGDQAFTE REAAEIMRDIGTAIQFLHSH NIAHRDVKPENLLYTSKEKD AVLKLTDFGFAKETTNALQ TPCYTPYYVAPEVLGPEKYD KSCDMWSLGVIMYILLCGF



Western blot analysis of lysates from wild type (WT) and MAPKAPK3 knockout (KO) HeLa cells, using [KO Validated] MAPKAPK3 Rabbit polyclonal antibody (STJ11100068) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081