

**Anti-RAB11FIP5 antibody (354-653) (STJ11100112)**  
STJ11100112

**GENERAL INFORMATION**




|                          |                    |
|--------------------------|--------------------|
| <b>Product Type</b>      | Primary antibodies |
| <b>Short Description</b> |                    |
| <b>Applications</b>      | WB/IF/ICC/ELISA    |
| <b>Host/Source</b>       | Rabbit             |
| <b>Reactivity</b>        | Human/Mouse/Rat    |

**PRODUCT PROPERTIES**

|                            |                                                                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                                              |
| <b>Clone ID</b>            |                                                                                                                                                         |
| <b>Concentration</b>       | Lot specific                                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                            |
| <b>Purification</b>        | Affinity purification                                                                                                                                   |
| <b>Dilution</b>            | WB:1:500-1:1000                                                                                                                                         |
| <b>Range</b>               | IF/CC:1:50-1:200<br>ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.                                                                                                        |
| <b>Isotype</b>             | IgG                                                                                                                                                     |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                          |

**TARGET INFORMATION**

|                           |                                                                                                                                                                                                                                                           |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 26056                                                                                                                                                                                                                                                     |
| <b>Gene Symbol</b>        | RAB11FIP5                                                                                                                                                                                                                                                 |
| <b>Uniprot ID</b>         | RFIP5_HUMAN                                                                                                                                                                                                                                               |
| <b>Immunogen</b>          |                                                                                                                                                                                                                                                           |
| <b>Immunogen Region</b>   | 354-653                                                                                                                                                                                                                                                   |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 354-653 of human RAB11FIP5 (NP_056285.1).                                                                                                                                   |
| <b>Immunogen Sequence</b> | HRSSISGSLPSSGSLQAVSS RFSEEGPRSTDDTWPRGSRSS NSSSEAVLGQEELSAQAKVL APGASHPGEEEGARLPEGKP VQVATPIVASSEAVAEKEGA RKEERKPRMGLFHHHHQGLS RSELGRRSSLGEKGGPILGA SPHHSSSGEEKAKSSWFGLR EAKDPTQKPSHPVKLSAA PVEGSPDRKQSRSSLSIALS SGLEKLKTVTSGSIQPVTAQ PQAGQMVDTKRLKDSAVLD |



Immunofluorescence analysis of A-431 cells using RAB11FIP5 rabbit polyclonal antibody (STJ11100112) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.

Immunofluorescence analysis of NIH/3T3 cells using RAB11FIP5 rabbit polyclonal antibody (STJ11100112) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081