

Anti-TBC1D21 antibody (1-336) (STJ11100183)

STJ11100183

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	161514
Gene Symbol	TBC1D21
Uniprot ID	TBC21_HUMAN
Immunogen	
Immunogen Region	1-336
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-336 of human TBC1D21 (NP_699187.1).
Immunogen Sequence	MTTSLPENSLSARQSASFIL VKRKPPIDKTEWDSFFDESG HLAKSRDFICVNILERGLHP FVRTEAWKFLTGYFSWQSSQ DERLTVDSMRKKNYKALCQM YEKIQPLENLHRNFTETRN NIARDIQKIYDKDPLGNVLI DKKRLEKILLLSYVCNTQAE YQQGFHEMMMLFQLMVEHDH ETFWLFQFFLQKTEHSCVIN IGVAKNLDMMLTLITFLDPV FAEHLKGKGAGAVQSLFPW