

Anti-CLTCL1 antibody (1600-1663) (STJ11100187)

STJ11100187

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8218
Gene Symbol	CLTCL1
Uniprot ID	CLH2_HUMAN
Immunogen	
Immunogen Region	1600-1663
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1600-1663 of human CLTCL1 (NP_009029.3).
Immunogen Sequence	IQVMREYLSKVDKLDALSL RKQEEHVTEPAPLVDFDFGH E