

Anti-ZHX2 antibody (488-837) (STJ11100192)

STJ11100192

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22882
Gene Symbol	ZHX2
Uniprot ID	ZHX2_HUMAN
Immunogen	
Immunogen Region	488-837
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 488-837 of human ZHX2 (NP_055758.1).
Immunogen Sequence	SDHRYRCQRGIVHITSESLA KDQLAIAASRHGRTHAYPD FAPQKFKEKTQGQVKILEDS FLKSSFPTQAEIDRLRVETK LSRREIDSWFSERRKLDSM EQAVLDSMGSGKKGQDVGAP NGALSRLDQLSGAQLTSSLP SPSPAIAKSQEQVHLLRSTF ARTQWPTPQEYDQLAAKTGL VRTEIVRWFKENRCLLKTGT VKWMEQYQHQPMADDHGYDA VARKATKPMESPKNNGDV



Western blot analysis of various lysates using ZHX2 Rabbit polyclonal antibody (STJ11100192) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081