

Anti-PRDX4 antibody (10-110) (STJ11100264)

STJ11100264

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	10549
Gene Symbol	PRDX4
Uniprot ID	PRDX4_HUMAN
Immunogen	
Immunogen Region	10-110
Specificity	A synthetic peptide corresponding to a sequence within amino acids 10-110 of human Peroxiredoxin 4 (PRDX4) (NP_006397.1).
Immunogen Sequence	TTPDHGRHRRLLPLLLFL LPAGAVQGWETEERPRTREE ECHFYGAGGQVYPGEASRVSV ADHSLHLSKAKISKAPPYWE GTAVIDGEFKEKLTIDYRGK Y




Immunofluorescence analysis of NIH/3T3 cells using Peroxiredoxin 4 (PRDX4) Rabbit polyclonal antibody (STJ11100264) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of Peroxiredoxin 4 (PRDX4) in paraffin-embedded mouse kidney using Peroxiredoxin 4 (PRDX4) Rabbit polyclonal antibody (STJ11100264) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of Peroxiredoxin 4 (PRDX4) in paraffin-embedded human liver cancer using Peroxiredoxin 4 (PRDX4) Rabbit polyclonal antibody (STJ11100264) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of Peroxiredoxin 4 (PRDX4) in paraffin-embedded rat liver using Peroxiredoxin 4 (PRDX4) Rabbit polyclonal antibody (STJ11100264) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081