

Anti-ADD3 antibody (400-500) (STJ11100292)

STJ11100292

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	120
Gene Symbol	ADD3
Uniprot ID	ADDG_HUMAN
Immunogen	
Immunogen Region	400-500
Specificity	A synthetic peptide corresponding to a sequence within amino acids 400-500 of human ADD3 (NP_001307523.1).
Immunogen Sequence	SGGTPIKIEDPNQFVPLNTN PNEVLEKRNKIREQNRYDLK TAGPQSLLAGIVVDKPPST MQFEDDDHGPPAPPNPFSHL TEGELEEKRTIERKQQGLE D