

Anti-PCIF1 antibody (1-200) (STJ11100302)

STJ11100302

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	63935
Gene Symbol	PCIF1
Uniprot ID	CAPAM_HUMAN
Immunogen	
Immunogen Region	1-200
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PCIF1 (NP_071387.1).
Immunogen Sequence	MANENHGSPREEASLLSHSP GTSNQSQCSPKPIRLVQDL PEELVHAGWEKCWSRRENRP YYFNRFTNQSLWEMPVLGQH DIVSDPLGLNATPLPQDSSL VETPPAENKPRKRQLSEEQP SGNGVKKPKIEIPVTPTGQS VPSSPSIPGTPTLKMWGTSP EDKQQAALLRPTEVYWDLDI QTNAVIKHRGPSEVLPPHPE



Western blot analysis of lysates from HeLa cells, using PCIF1 Rabbit polyclonal antibody (STJ11100302) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 20s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081