

Anti-PIK3R1/PIK3R2/PIK3R3 antibody (400 aa) (STJ11100308)

STJ11100308

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8503 5296 5295 PIK3R3 PIK3R2 P55G_HUMAN P85B_HUMAN P85A_HUMAN
Immunogen	400 aa
Immunogen Region	400 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 400 to 500 of human PIK3R1/PIK3R2/PIK3R3 (NP_852664.1).
Immunogen Sequence	



Western blot analysis of various lysates using PIK3R1/PIK3R2/PIK3R3 Rabbit pAb (STJ11100308) at 1:500 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (STJ5000858) at 1:10000 dilution.
Lysates/proteins: 25 Mu g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit
Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081