

Anti-PALM antibody (1-387) (STJ11100324)

STJ11100324

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5064
Gene Symbol	PALM
Uniprot ID	PALM_HUMAN
Immunogen	
Immunogen Region	1-387
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-387 of human PALM (NP_002570.2).
Immunogen Sequence	MEVLAAETTSQGERLQAIAE KRKRQAEIENKRRQLEDERR QLQHLKSKALRRERWLLGTP SSASEGDEDLRRQMQRDEQK TRLLEDSVSRLEKEIEVLER GDSAPATAKENAAAPSPVRA PAPSPAKEERKTEVWMNSQQ TPVGTPKDKRVSNTPLRTVD GSPMMKAAMYSVEITVEKDK VTGETRVLSSTTLPRQPLP LGIKVYEDETKVHVAVDGTAE ENGIHPLSSSEVDELIHKA