

Anti-RASA2 antibody (700-849) (STJ11100328)

STJ11100328

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5922
Gene Symbol	RASA2
Uniprot ID	RASA2_HUMAN
Immunogen	
Immunogen Region	700-849
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 700-849 of human RASA2 (NP_006497.2).
Immunogen Sequence	LCRVSRCNQNRLSFYHPSVY LNGNWLCQETGENTLGCKP CTAGVPADIQIDIDEDRETE RIYSLFTLSLLKLQKMEEAC GTIAVYQGPQKEPDDYSNFV IEDSVTTFKTIQIKSIIEK LDEPHEKYRKKRSSSAKYGS KENPIVGKAS