

Anti-TAF10 antibody (100 aa) (STJ11100340)

STJ11100340

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6881
Gene Symbol	TAF10
Uniprot ID	TAF10_HUMAN
Immunogen	
Immunogen Region	100 aa
Specificity	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human TAF10 (NP_006275.1).
Immunogen Sequence	PSAANGDVKPVVSSTPLVDF LMQLEDYTP TIPDAVTGYL NRAGFEASDPRIIRLISLAA QKFISDIANDALQHCKMKGT ASGSSRSKSKDRKYLTLMED LTPALSEYGINVKKPHYFT