

Anti-GCFC2 antibody (250-480) (STJ11100341)

STJ11100341

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6936
Gene Symbol	GCFC2
Uniprot ID	GCFC2_HUMAN
Immunogen	
Immunogen Region	250-480
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 250-480 of human C2orf3 (NP_003194.3).
Immunogen Sequence	GSSKVKKFDTSISFPPVNLE IIKKQLNTRLTLQETHRSH LREYEKYVQDVKSSKSTIQN LESSNQALNCKFYKSMKIY VENLIDCLNEIKINIEIES SMHALLLKQAMTFMKRRQDE LKHESTYLQQLSRKDETSTS GNFSVDEKTQWILEEIESRR TKRRQARVLSGNCNHQEGTS SDELPSAEMIDFQKSQGD I LQKQKKVFEEVQDDFCNIQN ILLKFQWREK