

Anti-PDE8B antibody (1-210) (STJ11100352)

STJ11100352

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8622
Gene Symbol	PDE8B
Uniprot ID	PDE8B_HUMAN
Immunogen	
Immunogen Region	1-210
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human PDE8B (NP_003710.1).
Immunogen Sequence	MGCAPSIHVSQSGVIYCRDS DESSSPROTTSVSQGPAAPL PGLFVQTDAAAIPPSRASG PPSVARVRRARTELGSGSSA GSAAPAAATSRGRRRHCCSS AEAETQTCYTSVKQVSSAEV RIGPMRLTQDPIQVLLIFAK EDSQSDGFWWACDRAGYRCN IARTPESALECFLDKHHEII VIDHRQTQNFDAEAVCRSIR ATNPSEHTVI