

Anti-IQCB1 antibody (500-600) (STJ11100362)

STJ11100362

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9657
Gene Symbol	IQCB1
Uniprot ID	IQCB1_HUMAN
Immunogen	
Immunogen Region	500-600
Specificity	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human IQCB1 (NP_001018864.2).
Immunogen Sequence	EERAQQHREALIAQISTNVE QLMKAPSLKEAEGKEPELFL SRSRPVAAKAKQAHLTTLKH IQAPWWKKLGEESGDEIDVP KDELSIELENLFIGGTKPP