

Anti-SEMA3E antibody (536-775) (STJ11100363)

STJ11100363

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9723
Gene Symbol	SEMA3E
Uniprot ID	SEM3E_HUMAN
Immunogen	
Immunogen Region	536-775
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 536-775 of human SEMA3E (NP_036563.1).
Immunogen Sequence	YCAWDGISCSRYPTGTHAK RFRFRQDVRHGNAAQQCFGQ QFVGDALDKTEEHLAYGIEN NSTLLECTPRSLQAKVWFV QKGRETRKEEVKTDDRNVKM DLGLLFLRLHKSDAGTYFCQ TVEHSFVHTVRKITLEWEE EKVEDMFNKDDEEDRHHMP CPAQSSISQGAKPWYKEFLQ LIGYSNFQRVEEYCEKVWCT DRKRKKLKMSPSKWKYANPQ EKKLRSKPEHYRLPRHTLD