

Anti-PHLDB1 antibody (650-950) (STJ11100382)

STJ11100382

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	23187
Gene Symbol	PHLDB1
Uniprot ID	PHLB1_HUMAN
Immunogen	
Immunogen Region	650-950
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 650-950 of human PHLDB1 (NP_001138231.1).
Immunogen Sequence	GRRPSRGLAGASGRSSEEPG VATQRLWESMERSDEENLKE ECSSTESTQQEHEDAPSTKL QGEVLALLEEERAQVLGHVEQ LKVRVKELEQQLESAREAE MERALLQGEREAERALLQKE QKAVDQLQEKLVALETGIQK ERDKEAELETETKLFEDLE FQQLERESRVEEERELAGQG LLRSKAELLRSIAKRKERLA ILDSQAGQIRAAVQESERL ARDKNASLQLQKEKEKLT