

Anti-ANKS1A antibody (700-950) (STJ11100389)

STJ11100389

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	23294
Gene Symbol	ANKS1A
Uniprot ID	ANS1A_HUMAN
Immunogen	
Immunogen Region	700-950
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 700-950 of human ANKS1A (NP_056060.2).
Immunogen Sequence	SVGEWLESIGLQQYESKLLL NGFDDVHFLGSNMEEQDLR DIGISDPQHRRKLLQAARSL PKVKALGYDGNSPPSVPSWL DSLGLQDYVHSFLSSGYSSI DTVKNLWELELVNVLKVQLL GHRKRIIASLADRPYEEPPQ KPPRFSQLRCQDLLSQTSPP LSQNDSCTEGRSADLLPPGD TGRRRHDSLHDAAPSRAER FRIQEEHREAKLTLRPPSLA APYAPVQSWQHQPKEKLIFE