

Anti-BRI3 antibody (1-100) (STJ11100396)

STJ11100396

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25798
Gene Symbol	BRI3
Uniprot ID	BRI3_HUMAN
Immunogen	
Immunogen Region	1-100
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human BRI3 (NP_056194.1).
Immunogen Sequence	MDHKPLLQERPPAYNLEAGQ GDYACGPHGYGAIPAAPPPP PYPYLVGTGIPTHPRVYNIH SRTVTRYPANSIVVGGCPV CRVGVLEDCFTFLGIFLAI