

Anti-FAM32A antibody (1-112) (STJ11100398)

STJ11100398

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26017
Gene Symbol	FAM32A
Uniprot ID	FA32A_HUMAN
Immunogen	
Immunogen Region	1-112
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-112 of human FAM32A (NP_054796.1).
Immunogen Sequence	MEAYEQVQKGPLKLGVAEL GVTKRKKKKDKDKAKLLEA MGTSKKNEEEKRRGLDKRTP AQAAFEKMQEKQMERILKK ASKTHKQRVEDFNRHLDLT EHYDIPKVSWTK