

Anti-SIPA1L1 antibody (1480-1700) (STJ11100399)

STJ11100399

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26037
Gene Symbol	SIPA1L1
Uniprot ID	S1L1_HUMAN
Immunogen	
Immunogen Region	1480-1700
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1480-1700 of human SIPA1L1 (NP_056371.1). TSKYLIGWKKPEGTINSVGF MDTRKRHQSDGNEIAHTRLR ASTRDLRASPKPTSKSTIEE DLKKLIDLESPTPESQKSFK FHALLSPQSPFPSTPTSRRA LHRTLSDESIYNSQREHFFT SRASLLDQALPNDVLFSTY PSLPKSLPLRRPSYTLGMKS LHGEFSASDSSLTIQETRR QPMPDPGLMPLPDTAADLDW SNLVDAAKAYEVQRASFFAA S
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081