

Anti-EPC2 antibody (528-807) (STJ11100401)

STJ11100401

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26122
Gene Symbol	EPC2
Uniprot ID	EPC2_HUMAN
Immunogen	
Immunogen Region	528-807
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 528-807 of human EPC2 (NP_056445.3).
Immunogen Sequence	FQPRLLNLQSDSEECTSRK PGQTVNNKRVSAASVALLNT SKNGISVTGGITEEQFQTHQ QQLVQMQRQQLAQLQQKQQS QHSSQQTHPKAQQGSSTSDCM SKTLDASAHFAASAVVSAP VPSRSEVAKEQNTGHNNING VVQPSGTSKTLYSTNMALSS SPGISAVQLVRTVGHTTTNH LIPALCTSSPQTLPMNNSCL TNAVHLNNVSVSPVNVHIN TRTSAPSPALKLATVAAS