

Anti-MRPS17 antibody (1-130) (STJ11100412)

STJ11100412

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51373
Gene Symbol	MRPS17
Uniprot ID	RT17_HUMAN
Immunogen	
Immunogen Region	1-130
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human MRPS17 (NP_057053.1).
Immunogen Sequence	MSVVRSSVHARWIVGKVIGT KMQKTAKVRVTRLVLDPYLL KYFNKRKTYFAHDAQQCTV GDIVLLRALPVPRAKHVKHE LAEIVFKVGVDPVTGKPC AGTTYLESPLSSETQLSKN LEELNISSAQ