

Anti-ZNF302 antibody (90-200) (STJ11100429)

STJ11100429

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55900
Gene Symbol	ZNF302
Uniprot ID	ZN302_HUMAN
Immunogen	
Immunogen Region	90-200
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 90-200 of human ZNF302 (NP_060913.2).
Immunogen Sequence	DEDSQPVTMEKVVKQSYEF SNSNKNLEYTECDTFRSTFH SKSTLSEPQNNSAEGNSHKY DILKKNLSSKSVIKSERING GKKLLNSNKGAAFNQSKSL TLPQTCNREKI