

Anti-AGXT2 antibody (42-300) (STJ11100440)

STJ11100440

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64902
Gene Symbol	AGXT2
Uniprot ID	AGT2_HUMAN
Immunogen	
Immunogen Region	42-300
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 42-300 of human AGXT2 (NP_001293102.1).
Immunogen Sequence	LHTKPRMPPCDFMPERYQSL GYNRVLEIHKEHLSPVVTAY FQKPLLLHQGHMEWLFDAEG SRYLDFFSGIVTSVGHCHP KVNAVAQKQLGRLWHTSTVF FHPPMHEYAEKLAALLPEPL KVIFLVNSGSEANELAMLMA RAHSNNIDIISFRGAYHGCS PYTLGLTNVGTYKMELPGGT GCQPTMCPDVFRGPWGGSHC RDSPVQTIRKCSAPDCCQA KDQYIEQFKDTLSTSVAKS