

Anti-ZCRB1 antibody (30-110) (STJ11100469)

STJ11100469

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	85437
Gene Symbol	ZCRB1
Uniprot ID	ZCRB1_HUMAN
Immunogen	
Immunogen Region	30-110
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 30-110 of human ZCRB1 (NP_149105.3).
Immunogen Sequence	FSKYGKVVKVTIMKDKDTRK SKGVAFILFLDKSAQNCTR AINNKKLFGRVIKASIAIDN GRAAEFIRRRNYFDKSKCYE C