

Anti-CHST14 antibody (180-280) (STJ11100478)

STJ11100478

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	113189
Gene Symbol	CHST14
Uniprot ID	CHSTE_HUMAN
Immunogen	
Immunogen Region	180-280
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 180-280 of human CHST14 (NP_569735.1).
Immunogen Sequence	LKMDHRSDLVFLADLRPEEI RYRLQHYFKFLVREPLERL LSAYRNKFGEIREYQQR YGA EIVRRYRAGAGPSPAGDDVT FPEFLRYLVDEDPERMNEHW M