

# Anti-NR1D1 antibody (280-360) (STJ11100546)

STJ11100546

## GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/IHC-P/IF/ICC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

## PRODUCT PROPERTIES


<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:500-1:2000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

## TARGET INFORMATION


<b>Gene ID</b>	9572
<b>Gene Symbol</b>	NR1D1
<b>Uniprot ID</b>	NR1D1_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	280-360
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 280-360 of human Rev-Erb Alpha/NR1D1 (NP_068370.1).
<b>Immunogen Sequence</b>	SPEPTVEDVISQVARAHREI FTYAHDKLGSSPGNFNANHA SGSPATTPHRWENQGCPPA PNDNNTLAAQRHNEALNGLR Q




Immunofluorescence analysis of PC-12 cells using Rev-Erb Alpha/NR1D1 Rabbit polyclonal antibody (STJ11100546) at dilution of 1:50 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using Rev-Erb Alpha/NR1D1 Rabbit polyclonal antibody (STJ11100546) at dilution of 1:50 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of Rev-Erb Alpha/NR1D1 in paraffin-embedded rat liver using Rev-Erb Alpha/NR1D1 Rabbit polyclonal antibody (STJ11100546) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of Rev-Erb Alpha/NR1D1 in paraffin-embedded mouse brain using Rev-Erb Alpha/NR1D1 Rabbit polyclonal antibody (STJ11100546) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081