

Anti-GTF2F1 antibody (390-517) (STJ11100599)

STJ11100599

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2962
Gene Symbol	GTF2F1
Uniprot ID	T2FA_HUMAN
Immunogen	
Immunogen Region	390-517
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 390-517 of GTF2F1 (NP_002087.2).
Immunogen Sequence	PSAEGGSTSTLRAAASKLE QGKRVSEMPAAKRLRLDTGP QSLSGKSTPQPPSGKTPNS GDVQVTEDAVRRYLTRKPMPT TKDLLKKFQTKKTGLSSEQT VNVLAQILKRLNPERKMIND KMHFSLKE



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081