

Anti-FKBP10 antibody (303-582) (STJ11100746)

STJ11100746

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	60681
Gene Symbol	FKBP10
Uniprot ID	FKB10_HUMAN
Immunogen	
Immunogen Region	303-582
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 303-582 of human FKBP10 (NP_068758.3).
Immunogen Sequence	FDSSYSRNHTYNTYIGQGYI IPGMDQGLQGACMGERRRIT IPPHLAYGENGTGDKIPGSA VLIFNVHVIDFHNPADVVEI RTLSRPSETCNETTKLGDFV RYHYNCSLLDGTQLFTSHDY GAPQEATLGANKVIEGLDTG LQGMCVGERRQLIVPPHLAH GESGARGVPGSAVLLFEVEL VSREDGLPTGYLFVWHKDP ANLFEDMDLNKDGEVPEEF STFIKAQVSEKGRLMPGQ



Western blot analysis of lysates from wild type (WT) and FKBP10 knockout (KO) 293T cells, using [KO Validated] FKBP10 Rabbit polyclonal antibody (STJ11100746) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1S.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081