

Anti-GSTK1 antibody (1-226) (STJ11100839)

STJ11100839

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	373156
Gene Symbol	GSTK1
Uniprot ID	GSTK1_HUMAN
Immunogen	
Immunogen Region	1-226
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-226 of human GSTK1 (NP_057001.1).
Immunogen Sequence	MGPLPRTVELFYDVLSPYSW LGFEILCRYQNIWNINLQLR PSLITGIMKDSGNKPPGLLP RKGLYMANDLKLLRHHLQIP IHFPKDFLSVMLEKGSLSAM RFLTAVNLEHPMLEKASRE LWMRVWSRNEDITEPQSILA AAEKAGMSAEQAQGLLEKIA TPKVKNQLKETTEAACRYGA FGLPITVAHVDGQTHMLFGS DRMELLAHLLGEKWMGPIPP AVNARL



Western blot analysis of lysates from wild type (WT) and GSTK1 knockout (KO) 293T cells, using [KO Validated] GSTK1 Rabbit polyclonal antibody (STJ11100839) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) [STJS000859] at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081