

Anti-LAMTOR1 antibody (1-161) (STJ11100859)

STJ11100859

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IF/CC:1:50-1:100 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55004
Gene Symbol	LAMTOR1
Uniprot ID	LTOR1_HUMAN
Immunogen	
Immunogen Region	1-161
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-161 of human LAMTOR1 (NP_060377.1).
Immunogen Sequence	MGCCYSSENESSDQDREERK LLLDPSSPPTKALNGAEPNY HSLPSARTDEQALLSSILAK TASNIIDVSAADSQGMEQHE YMDRARQYSTRLAVLSSSLT HWKKLPPLSLTSQPHQVLA SEPIPFSDLQQVSRIAAYAY SALSQIRVDAKEELVWQFGI P



Western blot analysis of lysates from wild type (WT) and LAMTOR1 knockout (KO) HeLa cells, using [KO Validated] LAMTOR1 Rabbit polyclonal antibody (STJ11100859) at 1:1000 dilution. Secondary antibody: HRP-Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081