

Anti-SKI antibody (700-800) (STJ11100892)

STJ11100892

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6497
Gene Symbol	SKI
Uniprot ID	SKI_HUMAN
Immunogen	
Immunogen Region	700-800
Specificity	A synthetic peptide corresponding to a sequence within amino acids 650-728 of human SKI (NP_003027.1).
Immunogen Sequence	CDKGCEAGRLRAKYSAQIED LQVKLQHAEDREQLRADLL REREAREHLEKVVKELQEQL WPRARPEAAGSEGAEELEP



Western blot analysis of extracts of mouse fetal muscle, using SKI antibody (STJ11100892). Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJUS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081