

Anti-NSUN5 antibody (127-429) (STJ11100922)

STJ11100922

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:1000
Range	IHC-P: 1:50-1:200
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	55695
Gene Symbol	NSUN5
Uniprot ID	NSUN5_HUMAN
Immunogen	
Immunogen Region	127-429
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 127-429 of human NSUN5 (NP_060514.1).
Immunogen Sequence	RPGPASQLPRFVRVNTLKTC SDDVVDYFKRQGFYQGRAS SLDDLRLKLGKHFLDPLMP ELLVFPAQTDLHEHPLYRAG HLILQDRASCLPAMLLDPPP GSHVIDACAAPGNKTSHLAA LLKNQGGKIFAFDLDAKRLAS MATLLARAGVSCCELAEDF LAVSPSPDPHYEHVYILLDP SCSSGSGMPSPRQLEEPGAGTP SPVRLHALAGFQQRALCHAL TFPSLQRLVYSTCSLCQEE




Immunohistochemistry analysis of paraffin-embedded Rat brain using NSUN5 antibody (STJ11100922) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded Mouse spinal cord using NSUN5 antibody (STJ11100922) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded Mouse kidney using NSUN5 antibody (STJ11100922) at dilution of 1:100 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of extracts of various cell lines, using NSUN5 antibody (STJ11100922) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081