

Anti-MAS1 antibody (246-325) (STJ11100947)

STJ11100947

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4142
Gene Symbol	MAS1
Uniprot ID	MAS_HUMAN
Immunogen	
Immunogen Region	246-325
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 246-325 of human MAS1 (NP_002368.1).
Immunogen Sequence	LLYLLYYEYWSTFGNLHHIS LLFSTINSSANPFIYFFVGS SKKKRFKESLKVLTAFKD EMQPRRQKDNCTVTETVV