

## Anti-Phospho-HNF4A-S304 antibody (STJ11100975)

STJ11100975

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	WB/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	3172
<b>Gene Symbol</b>	HNF4A
<b>Uniprot ID</b>	HNF4A_HUMAN
<b>Immunogen</b>	SDPGK
<b>Immunogen Region</b>	
<b>Specificity</b>	A synthetic phosphorylated peptide around S304 of human HNF4A (NP_849180.1).
<b>Immunogen Sequence</b>	SDPGK