

Anti-Phospho-DDX39A-Y138 antibody (STJ11101054)

STJ11101054

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10212
Gene Symbol	DDX39A
Uniprot ID	DX39A_HUMAN
Immunogen	SKYMP
Immunogen Region	
Specificity	A synthetic phosphorylated peptide around Y138 of human DDX39A (NP_005795.2).
Immunogen Sequence	SKYMP